

PROTOCO₂L TOUCH[®] Colon Insufflator

INDICATIONS AND USAGE:

The PROTOCO₂L TOUCH[®] COLON INSUFFLATOR administers and regulates carbon dioxide as a distention media to the colon during CT Colonography (CTC or Virtual Colonoscopy) and conventional Colonoscopy.

IMPORTANT SAFETY INFORMATION / CONTRAINDICATIONS FOR USE:

The PROTOCO₂L TOUCH COLON INSUFFLATOR should be used only when colon insufflation is indicated, and should therefore not be used for any other treatments. It should only be used under the direct guidance of a physician experienced in colon insufflation.

This device is **contraindicated** for hysteroscopic insufflation, i.e., it must not be used for intrauterine distention.

This product should not be used in patients with known or suspected colonic perforation or toxic megacolon. It should not be used within 6 days of large forceps or “hot” biopsy, or snare polypectomy.

Do not use this product in a colostomy stoma.

Do not use this product following recent rectal surgery or low rectal anastomosis, or when proctitis or other rectal conditions such as inflammatory or neoplastic diseases are suspected.

The PROTOCO₂L TOUCH Operator’s Manual may be obtained by contacting Bracco Diagnostics Inc. Professional Services Department at 800-257-5181, Option 2.

PROTOCO₂L TOUCH is manufactured for Bracco Diagnostics Inc., Princeton, NJ 08540.

PROTOCO₂L TOUCH is a registered trademark of E-Z-EM, Inc.

Bracco Diagnostics Inc.
510 Carnegie Center
Suite 300
Princeton, NJ 08540 USA
Phone: 609-514-2200
Toll-Free: 1-800-631-5245 (U.S. only)
Fax: 609-514-2424

© 2025 Bracco Diagnostics Inc. All Rights Reserved.