



MR Suite
SOLUTIONS



I/Dx Suite



Nuclear Suite



CT Suite

Moving to the **rhythm** of your workflow

An integrated system designed to work with
you to achieve optimal image quality



LIFE FROM INSIDE

EmpowerMR[®]

Contrast Injection System

YOU



EmpowerMR[®] eliminates RF transients



A state-of-the art contrast injection system compatible with scanners up to 7T

The EmpowerMR[®] contrast injection system offers a full range of unique features to optimize safety while helping technologists achieve optimal images

ENHANCING YOUR PRODUCTIVITY

Unique features	<i>to increase efficiency</i>
Intuitive, user-friendly interface	<i>Easy operation</i>
Compact, rotating head	<i>More operational flexibility</i>
Dual initialize and replace syringe	<i>Increases productivity</i>
Non-magnetic injector head construction	<i>No limits on proximity to scanner</i>

ADVANCING PATIENT SAFETY

Unique features	<i>to protect your patients</i>
Battery-free hydraulic drive system	<i>Continuous power supply avoids risk of battery failure</i>
Tilt lockout	<i>Helps to minimize risk of air embolism</i>
Color-coded, illuminating handle	<i>Provides notification that injector ram is moving</i>
Overpressure during KVO	<i>Pop-up message on remote notifies of possible occlusions</i>
GFR calculator	<i>Helps identify patients with compromised renal function</i>

IRiSMR[®] DATA NETWORKING APPLICATION PROMOTES INFORMED DECISION-MAKING

- Captures and consolidates essential contrast and syringe use data
- Generates management reports to help determine budgets, control cost, and enhance workflow
- Automatically archives injection information
- Exports data to Microsoft[®] Office Suite for analysis and review*
- Simplified summary screen provides real-time data for timely decision-making

* Network version

EmpowerMR[®]
Contrast Injection System

Moving to the **rhythm** of your workflow

EmpowerMR[®]
Contrast Injection System

ORDERING INFORMATION

Basic System Description

Catalog No.



EmpowerMR[®] Contrast Injection System

Includes: Injector head with hydraulic controller, touch screen remote (Windows[®] operating system based) with desktop mount, electro-hydraulic interconnect, communication and power cables (includes pedestal mount for injector head)


102901



Fast-Load[™] Syringes

2 disposable 100 mL Fast-Load Syringes with 96" (y) low-pressure connecting tube and bottle spikes

100301



IRiSMR[®] Data Networking Application

IRiSMR[®] Remote Viewer

100401

100404 (no charge)

BRACCO SETS THE PACE WITH INNOVATIVE **SOLUTIONS** TO ADDRESS YOUR IMAGING CHALLENGES

The EmpowerMR Injector System is used to inject MR contrast media followed by a saline flushing solution into the human vascular system in conjunction with magnetic imaging resonance (MRI) procedures.

Bracco Diagnostics Inc.
Distributed by Bracco Diagnostics Inc.
107 College Road East
Princeton, New Jersey 08540
Phone: (800) 631-5245
Fax: (609) 514-2424

Customer Service:
1-877-BRACCO 9 (1-877-272-2269)

Professional Services:
1-800-257-5181 (option 1)

Email: BRACCO@diag.bracco.com
Website: usa.braccoimaging.com

ACIST, EmpowerMR and IRiSMR are registered trademarks of ACIST Medical Systems, Inc.
Fast-Load is a trademark of ACIST Medical Systems, Inc.
Not all products available in all global markets.

ACIST Medical Systems, Inc.
Manufactured by ACIST Medical Systems, Inc.
Corporate Headquarters
7905 Fuller Road
Eden Prairie, Minnesota 55344
Phone: (952) 941-3507
Fax: (952) 942-4648

Toll-free in U.S.: 1-888-667-6648
Website: www.acist.com

© 2010 Bracco Diagnostics Inc. All Rights Reserved.
Printed in USA 10-062210D



Committed to Science,
Committed to You[™].

